

Emergency First Aid box with Eyewash Bottles

Where water is not immediately available the BROEN-LAB Emergency First Aid box can provide assistance quickly. Its contents mean treatment can begin immediately saving time required to move to a fixed emergency shower, before medical treatment.

The BROEN-LAB Emergency First Aid box can have its contents configured to meet local needs making it both versatile and relevant to the hazards present:

- Rinsing bottle with eyewash solution (to rinse out foreign bodies and for continuous rinsing until medical treatment)
- Eyewash bottle with pH neutral solution (for rinsing in the case of injury due to acids or alkaline solutions)
- The rinsing bottles are designed to rinse both eyes simultaneously maximizing and speeding the removal of particles and chemicals
- Wound cleaning wipes
- Wound dressings
- When removed, self-opening plaster dispensers
- Cooling gel



Emergency First Aid box with Eyewash Bottles

Details:

Robust plastic box with clear perspex front.

Emergency plan

No.: P7 33 005 10 00 00 2

Dimensions:

253 x 430 x 92 mm (WxHxD)

Standards:

Conforms to EN 15154 Part 4

Content:

- Safety symbol set with user instructions and test label
No. 17 001.970
- Sterile eyewash solution 0.5 l
No.: P7 17 002 10 00 00 0
- pH neutral eyewash solution 0.5 l
No.: P7 17 003 10 00 000
- Sterile Quick-Stop wound dressing
No.: P7 53 008 10 00 00 0
- QuickClean wound cleaner
No.: P7 53 007 10 00 00 0
- QuickCool burn gel
No.: P7 53 005 10 00 00 0
- 2 x QuickFix plasters: elastic
No.: P7 53 006 10 00 00 0

Optional Plasters:

- QuickFix, Waterproof
No.: P7 53 005 10 00 00 0
- QuickFix, Detectable
No.: P7 53 004 10 00 00 0
- QuickFix, Elastic, long
No.: P7 53 009 10 00 00 0

The content of the box must be checked regularly to ensure that it is complete and that expiration dates have not passed.

